

UHF Firmware Upgrade Guide

Software Release 13.02

Zurich
Instruments

Software Release 13.02 for the UHF Lock-in Amplifier

Zurich Instruments software release 13.02 consists of both software that runs on your PC and software (firmware) that runs on UHF Lock-in Amplifier itself. In order for the PC software in Release 13.02 to operate correctly it is also necessary to update the UHF's firmware. This step only has to be performed once. The aim of this guide is to explain how to upgrade the UHF's firmware.

Preparation

In order to upgrade the UHF's firmware, you must first take the following steps:

1. [Download](#) and install the appropriate version (32bit/64bit) of "LabOne 13.02" on your PC. Administrator rights are necessary for the installation. Please see the UHF User Manual.
2. Start the UHF Lock-in Amplifier.
3. Connect the UHF to the PC with the LabOne 13.02 installation via USB cable.

Starting the UHF Firmware Upgrade Utility

The UHF Firmware Upgrade Utility is the program used to perform a UHF firmware upgrade, it's a GUI (Graphical User Interface) included in a standard 13.02 LabOne installation.




Figure 1: Starting the Firmware Upgrade Utility via the Windows Start Menu.

In order to start the Firmware Upgrade Utility:

1. Click the Windows Icon in order to open the Start Menu
2. Click on the entry “All Programs” and navigate to and click on the “Zurich Instruments” folder
3. Under the “Servers” folder,click on the entry “Firmware upgrade”. This launches the Firmware Upgrade Utility.

Note, it’s not necessary to have administrator rights in order to start or use the UHF Firmware Upgrade Utility.

Using the UHF Firmware Upgrade Utility

Important: Do not disconnect the USB cable to the UHF or power-cycle the UHF whilst performing any of the following steps.

Upon starting the Firmware Upgrade Utility it should detect the device that is connected to the PC via USB. The device ID is displayed next to “Device:”.




Figure 2: The UHF Firmware Upgrade Utility upon start-up.

Select the device you would like to upgrade

Select which device you would like to upgrade via the pull-down menu. If no device is listed, please try the following steps:

1. Ensure that the USB cable is properly connected.
2. Try power-cycling the device.
3. Click the “Refresh” button.

Verify the firmware of the connected device

Click the “Verify” button to check the version of the current firmware.




Figure 3: Verifying the UHF firmware version.

Program the new firmware to the connected device

Click the “Program” button to install the UHF’s firmware.

Important: After verifying the firmware or programming the new firmware it is necessary to power-cycle the UHF for it to resume normal operation.




Figure 4: It is necessary to power-cycle the device after verifying or programming the firmware.

Close the UHF Firmware Upgrade Utility

Click the “Exit” button to close the UHF Firmware Upgrade Utility.

Support: If you encounter any issues whilst upgrading the UHF firmware, please contact: support@zhinst.com.